



# Lynn Pierri DDS, MS

Caring Without Compromise™

## SCAN STENT FABRICATION FOR MATERIALISE SIMPLANT

**Full Arch:** Holes drilled in the anterior section of the arch are placed midway between the incisal edge and cingulum. Holes drilled in the posterior teeth are placed in the center of the tooth.



**Posterior Tooth-Supported Stent:** Holes drilled in the posterior teeth are placed in the center of the tooth: First hole is 4mm from the tooth, then, 7mm for subsequent implants for biscuspidisation, or 8mm for molars.



**Anterior Tooth-Supported Stent:**  
Holes drilled in the anterior section of the arch are placed midway between the incisal edge and cingulum.

